

→ Autoclaved — To remove other microbes.

15 atm, 121°C , 15 min

⊗ Temperature, pH is needed to grow bacteria too.

Bio Lab →

8th January 2024

LB Media →

100ml water


Tryptone — 1 gm

Yeast extract — 0.5 gm


NaCl — 1 gm

Growth curve →

We would grow E. coli




If we are given a sample, how can we estimate the number of bacteria in sample?



Serial dilution

10x each time



→ Spread on petri

→ Count CFUs

~~If you~~ If you have > 300 , difficult to count.

If you have < 30 , too little (statistically bad)

~~$50 \rightarrow 10^{-6} \text{ L}$~~

~~$1 \rightarrow$~~


~~$(4 \text{ L} = 50 \times 10^6) 10^{-8}$~~

~~$1 \text{ L} = 50 \times 10^9$~~

Now say for about $10^3/\text{L}$ num of bacteria?

The number is too small for the previous techniques.

→ We centrifuge, take pellet and dissolve it in smaller volume, then repeat.



→ We could also use molecular ~~se~~ sieve paper to filter solution, dissolve what is on paper in lower volume, then repeat prev method.

If we want quicker results, we take OD.




$$\lambda = 600 \text{ nm}$$

OD 1 $\equiv 10^9 / \text{mL}$ \hookrightarrow Culture

We O.D the culture at fixed intervals, and use the fixed (determined solutions) to calculate.

Today's ~~exper~~ experiment \rightarrow

- (*) Prepare LB agar
- (*) Autoclave
- (*) Plug.



Generation time - 20 - 25 min (*E. coli*)

Bigger flask \rightarrow 50mL media

Smaller flask \rightarrow 20mL media